

AMPri Handelsgesellschaft mbH Benzstraße 16 21423 Winsen (Luhe) - Germany Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: Description: 05561-Y-XXL MED-COMFORT non-woven gown with PE-coating



SIZES	
	XXL
width	150 cm
length	140 cm

PRODUCT DES	RODUCT DESCRIPTION													
material	Ø	Polypro- pylene (PP)	V	Polyethy- lene (PE)		SMS								
colour		blue [B]		green [G]		orange [O]		red [R]		white [W]	V	yellow [Y]		light-blue [LB]
		dark-blue [DB]		dark-green [DG]										
closure	V	to tie		zipper		push bottons		velcro						
characteristics	V	tricot-cuff		elastic wrist- cuff	V	latex-free		hood		stand-up collar	N	partially laminated		sewed PE- oversleeves
		shortenable by perforation		blocked		dust-cover		embossed		short sleeves		with cardboard	V	single folded

REGULATORY AFFAIRS													
PPE-Regulation (EU) 2016/425		category I		category II		category III	V	no PPE- article					CE
MD-Regulation (EU) 2017/745	V	class I		class II		class III		sterile		measuring fundtion		no medical device	
Food Contact (EG) 1935/2004		acidic foods		aqueous foods		fatty foods		alcoholic foods		dry foods	V	not approved for food-contact	

STANDARDS													
tested in	V	DIN EN	V	ISO 16603	V	ISO 16604	V	ISO 22610	V	ISO 22611	V	ISO 22612	
accordance to		14126											
		EN ISO		EN 13034									
	Ŀ	13688											
	Ŀ												



AMPri Handelsgesellschaft mbH Benzstraße 16 21423 Winsen (Luhe) - Germany Telefon: +49 (0)4171 / 8480-0 Homepage: www.ampri.de e-mail: info@ampri.de

Technical Data Sheet

Article-No.: Description:

No.: 05561-Y-XXL

: MED-COMFORT non-woven gown with PE-coating



LOGISTIC DATA SUBPACKING							
general information							
material		bag					
pieces per subpacking	5	10					
GTIN subpacking		4044941022512					
measures & sizes							
length		430 mm					
width		330 mm					
height	-						
weights							
size	net weight	brut weight					
XXL 830 g 835 g							

LOGISTIC DATA PALETTE							
general information							
kind of palette euro-palette							
measures & sizes							
cartons per layer 4							
layer per palette		6					
palette height	palette height 201						
weights							
size	brut weight						
XXL 212 kg 237 kg							

LOGISTIC DATA OUTER PACKING

general information							
material		carton					
subpackings per outer	packing	10					
GTIN outer packing		4044941022529					
measures & sizes							
length		530 mm					
width		280 mm					
height		310 mm					
weights							
size	net weight	brut weight					
XXL	8.350 g	8.850 g					



Technical Data Sheet

Article-No.: 05561-Y-XXL

Description: MED-COMFORT non-woven gown with PE-coating



WARNINGS AND SAFETY INFORMATION

storage / expiry date	Store product in original packaging in a cool and dry place without additional weight, protect from direct sunlight. Do not store near ozone sources (laser printers, copiers). The actual expiry time in use cannot be specified in general terms, as it depends on the general conditions of use. An individual risk assessment must be carried out in each case. The expiry date - valid for proper storage - is stated on the packaging.
use and control	Only ever use the product for its intended use. A check/risk assessment must be carried out to ensure that the product is suitable for the intended use. Ensure that the size matches that of the user. Do not make any changes to the product. Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.). This product is intended for single use only
disposal	The used product must be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused products can be disposed of with household waste.
warnings	Prolonged wearing can lead to heat stress.

	This product is flammable - keep away from fire.
donning and doffin	g WAY OF PUTTING ON:
instructions	- Make sure that the size matches that of the user. Do not make any changes to the product.

domining and doming	WAT OF POTTING ON:
instructions	- Make sure that the size matches that of the user. Do not make any changes to the product.
	- Make sure that the product has no defects and is in good condition (no holes, unstitched parts, etc.).
	- Open the gown and put it on. Make sure that the garment does not tear. Close the gown by wrapping the waistband
	around your body and fastening it with a bow.
	REMOVAL METHOD:
	- Undo the waist and collar ties and remove the gown.

rev-no.: 8 date: 04.04.2025

changes and errors excepted