

Technical Data Sheet

Article-No.: **701190**
Description: **MED-COMFORT SENSITIVE**
shampoo cap
impregnated with quality lotion with
conditioner



DIMENSIONS				
length	width	height	diameter	capacity
-	-	-	-	-

PRODUCT DESCRIPTION							
material	☐ Polypropylene (PP)	☐ Polyethylene (PE)	☐ PP non woven	☐ PP non woven with PE coating	☐ paper	☒ 50% viscose + 50% polyester with PE coating	
colour	☐ blue [B]	☐ cedro [C]	☐ dark blue [DB]	☐ green [G]	☐ green opac [GO]	☐ purple [L]	☐ lightblue [LB]
	☐ Lime [LIME]	☐ Mix [MIX] black, cedro, purple, pink, lightblue	☐ orange [O]	☐ pink [P]	☐ red [R]	☐ rose [ROSE]	☐ black [SCH]
	☐ violett [V]	☒ white [W]	☐ wine red [WR]	☐ yellow [Y]	☐ silver	☐ natural	☐ bordeaux
characteristics	☒ single use	☐ single bladed	☐ 1 ply	☐ lint free	☐ with white lid	capacity: basin: approx. 18 litres bag: approx. 9 litres	
ingredients	☐ 50% Urea (CAS 57-13-6)	☐ 50% water	☐ 70 % isopropyl alcohol				
remark	Light and nourishing hair wash, without water. Alcohol-free, dermatologically tested. Mild and cleansing lotion.						

REGULATORY AFFAIRS							
PPE-Regulation (EU) 2016/425	☐ Category I	☐ Category II	☐ Category III	☒ no PPE-article			
#BEZUG!	☐ Class I	☐ Class II	☐ Class III	☐ sterile	☐ measuring function	☒ no medical device	
Food Contact (EG) 1935/2004	☐ acidic foods	☐ aqueous foods	☐ fatty foods	☐ alcoholic foods	☐ dry foods	☒ not approved for food-contact	

Technical Data Sheet

Article-No.: **701190**

Description: **MED-COMFORT SENSITIVE**
shampoo cap
impregnated with quality lotion with
conditioner



LOGISTIC DATA SUBPACKING		
general information		
material	bag	
pieces per subpacking	1	
GTIN subpacking blue [B]	-	
GTIN subpacking cedro [C]	-	
GTIN subpacking dark blue [DB]	-	
GTIN subpacking green [G]	-	
GTIN subpacking green opac [GO]	-	
GTIN subpacking purple [L]	-	
GTIN subpacking lightblue [LB]	-	
GTIN subpacking Lime [LIME]	-	
GTIN subpacking Mix [MIX]	-	
GTIN subpacking orange [O]	-	
GTIN subpacking pink [P]	-	
GTIN subpacking rot [R]	-	
GTIN subpacking rose [ROSE]	-	
GTIN subpacking black [SCH]	-	
GTIN subpacking violet [V]	-	
GTIN subpacking white [W]	4044941041551	
GTIN subpacking wine red [WR]	-	
GTIN subpacking yellow [Y]	-	
GTIN subpacking silver	-	
GTIN subpacking natur	-	
GTIN subpacking	-	
measures & size		
length	220 mm	
width	130 mm	
height	-	
weights		
size	net weight	brut weight
uni	130 g	135 g

LOGISTIC DATA OUTER PACKING		
general information		
material	carton	
subpackings per outer packing	24	
GTIN oute packing blue [B]	-	
GTIN oute packing cedro [C]	-	
GTIN oute packing dark blue [DB]	-	
GTIN oute packing green [G]	-	
GTIN oute packing green opac [GO]	-	
GTIN oute packing purple [L]	-	
GTIN oute packing lightblue [LB]	-	
GTIN oute packing Lime [LIME]	-	
GTIN oute packing Mix [MIX]	-	
GTIN oute packing orange [O]	-	
GTIN oute packing pink [P]	-	
GTIN oute packing rot [R]	-	
GTIN oute packing rose [ROSE]	-	
GTIN oute packing black [SCH]	-	
GTIN oute packing violet [V]	-	
GTIN oute packing white [W]	4044941041568	
GTIN oute packing wine red [WR]	-	
GTIN oute packing yellow [Y]	-	
GTIN oute packing silver	-	
GTIN oute packing natur	-	
GTIN oute packing	-	
measures & size		
length	220 mm	
width	200 mm	
heigh	335 mm	
weights		
size	net weight	brut weight
uni	3.240 g	3.740 g

LOGISTIC DATA PALLETTE		
general information		
kind of palette	euro-palette	
measures & size		
cartons per layer	11	
layer per palette	7	
height of the palette	250 cm	
weights		
size	net weight	brut weight
uni	288 kg	313 kg

Technical Data Sheet

Article-No.: **701190**

Description: **MED-COMFORT SENSITIVE**
shampoo cap
impregnated with quality lotion with
conditioner



WARNINGS AND SAFETY INFORMATION

storage	Store items in original packaging in a cool and dry place without additional weight, protect from direct sunlight.
use and control	Massage the cap onto the head for 2 to 3 minutes. Remove the cap, dry the hair and shape it. Rinsing is not necessary. For more convenience, the pack with cap can be heated in the microwave (30 sec, 600 W, always check the temperature). If the product is heated in the microwave before use, it is necessary to check the temperature after heating and before use.
cleaning instructions	disposable - cleaning is not intended for this article.
disposal	The product must be disposed in accordance with the disposal regulations of the local waste disposal company. Unused products can be disposed of with household waste. Do not dispose of in the toilet.
warnings / safety information	For external use only. Avoid contact with the eyes.

rev-no.: 2
date: 07. 01.2025

changes and errors excepted