

Technical Data Sheet

Article-No.: **01168**
Description: **EPIDERM PROTECT Purple**
Nitrile examination glove
purple, non sterile, powder free



| THICKNESS | | |
|-----------|---------|------------|
| cuff | palm | fingertips |
| 0,05 mm | 0,06 mm | 0,07 mm |

| PRODUCT DESCRIPTION | | | | | | | |
|---------------------|---|---|---------------------------------------|---|--|---|--|
| material | <input type="checkbox"/> Latex | <input checked="" type="checkbox"/> Nitrile | <input type="checkbox"/> Vinyl | <input type="checkbox"/> Vinyl-Nitrile-mixture | <input type="checkbox"/> Polyethylen (PE) | <input type="checkbox"/> TPE | <input type="checkbox"/> cotton |
| colour | <input type="checkbox"/> white | <input type="checkbox"/> blue | <input type="checkbox"/> black | <input type="checkbox"/> mint | <input checked="" type="checkbox"/> purple | <input type="checkbox"/> mix | <input type="checkbox"/> metallic blue |
| characteristics | <input type="checkbox"/> pre-powdered | <input checked="" type="checkbox"/> powder free | <input type="checkbox"/> sterile | <input checked="" type="checkbox"/> non sterile | <input checked="" type="checkbox"/> ambidextrous | <input type="checkbox"/> fits hand-specific | <input checked="" type="checkbox"/> accelerator-free |
| surface | <input type="checkbox"/> fully textured | <input checked="" type="checkbox"/> finger textured | <input type="checkbox"/> non textured | <input type="checkbox"/> embossed | <input checked="" type="checkbox"/> chlorinated inside | | |

| SIZES | | | | | | | |
|--------|----------|------------|------------|-------------|-------------|-------------|--------------|
| | XS (5-6) | S (6-7) | M (7-8) | L (8-9) | XL (9-10) | XXL (10-11) | XXXL (11-12) |
| width | ≤ 80 mm | 80 ± 10 mm | 95 ± 10 mm | 110 ± 10 mm | 115 ± 10 mm | - | - |
| length | ≥ 240 mm | ≥ 240 mm | ≥ 240 mm | ≥ 240 mm | ≥ 240 mm | - | - |

| REGULATORY AFFAIRS | | | | | | | |
|------------------------------|--|---|--|---|---|--|----|
| PPE-Regulation (EU) 2016/425 | <input type="checkbox"/> Category I | <input type="checkbox"/> Category II | <input checked="" type="checkbox"/> Category III | <input type="checkbox"/> no PPE-article | | | CE |
| MD-Regulation (EU) 2017/745 | <input checked="" type="checkbox"/> Class I | <input type="checkbox"/> Class II | <input type="checkbox"/> Class III | <input type="checkbox"/> sterile | <input type="checkbox"/> measuring function | <input type="checkbox"/> no medical device | |
| Food Contact (EG) 1935/2004 | <input checked="" type="checkbox"/> acidic foods | <input checked="" type="checkbox"/> aqueous foods | <input checked="" type="checkbox"/> fatty foods | <input checked="" type="checkbox"/> alcoholic foods | <input checked="" type="checkbox"/> dry foods | <input type="checkbox"/> not approved for food-contact | |

| STANDARDISATION | | | | | | |
|------------------------------|---------------------|---------------------------------|-----------------|---------------------|-------------------------------|-------------|
| EN 388:2016 Mechanical Risks | abrasion resistance | blade cut resistance Coupe-Test | tear resistance | puncture resistance | blade cut resistance TDM-Test | impact test |
| Level | not applicable | | | | | |

| | | | | | | |
|---------------------------------------|--|-------------|-------|-----------------|-------------|-----------------------|
| EN 374-1:2016 Chemical Risks | chemical | code letter | level | permeation time | degradation | ISO 374-1/Type B |
| | Sodium hydroxide 40% | K | 6 | > 480 min | 2,8 | |
| | Hydrogen Peroxide 30% | P | 4 | > 120 min | 32,1 | |
| EN 374-4:2013 Degradation | Formaldehyde 37% | T | 3 | > 60 min | 24,5 | KPT |
| | | | | | | |
| EN 374-5:2016 microorganism tightness | The glove is tight against microorganisms (viral, bacteria and fungi). Test according to ISO 16604 - method B. | | | | | EN ISO 374-5:2016 |

| | | | | | | |
|-------------------------------------|---|--|--|--|--|--|
| EN ISO 21420:2020 protective gloves | The glove meets the requirements according to EN ISO 21420:2020 | | | | | |
|-------------------------------------|---|--|--|--|--|--|

| | | | | | | |
|--------------------------------------|--|--|--|--|--|---------|
| EN 455 medical gloves for single use | The glove meets the requirements according to EN 455-1, EN 455-2, EN 455-3, EN 455-4 | | | | | EN 455 |
| EN 455-1:2000 freedom from holes | The glove has an AQL < 1.5 in regards to the water retention test (sampling inspection in acc. to ISO 2859-1, general Inspection Level I) | | | | | AQL 1.5 |

| | | | | | | |
|-----------------------------------|----------------|--|--|--|--|--|
| EN 16350 electrostatic properties | not applicable | | | | | |
|-----------------------------------|----------------|--|--|--|--|--|

Technical Data Sheet

Article-No.: **01168**
Description: **EPIDERM PROTECT Purple**
Nitrile examination glove
purple, non sterile, powder free



LOGISTIC DATA SUBPACKING

| generell information | | |
|--------------------------|---------------|--------------|
| material | carton | |
| pieces per subpacking | 100 | |
| EAN subpacking size XS | 4044941723624 | |
| EAN subpacking size S | 4044941723648 | |
| EAN subpacking size M | 4044941723662 | |
| EAN subpacking size L | 4044941723686 | |
| EAN subpacking size XL | 4044941723709 | |
| EAN subpacking size XXL | - | |
| EAN subpacking size XXXL | - | |
| PZN subpacking size XS | 17896532 | |
| PZN subpacking size S | 17896503 | |
| PZN subpacking size M | 17896495 | |
| PZN subpacking size L | 17896489 | |
| PZN subpacking size XL | 17896526 | |
| PZN subpacking size XXL | - | |
| PZN subpacking size XXXL | - | |
| measures & size | | |
| length | 240 mm | |
| width | 120 mm | |
| height | 50 mm | |
| weights | | |
| size | net weight | gross weight |
| XS | 310 g | 370 g |
| S | 340 g | 400 g |
| M | 370 g | 430 g |
| L | 400 g | 460 g |
| XL | 430 g | 490 g |
| XXL | - | - |
| XXXL | - | - |

LOGISTIC DATA OUTER PACKING

| generell information | | |
|-------------------------------|---------------|--------------|
| material | carton | |
| subpackings per outer packing | 10 | |
| EAN outer packing size XS | 4044941723631 | |
| EAN outer packing size S | 4044941723655 | |
| EAN outer packing size M | 4044941723679 | |
| EAN outer packing size L | 4044941723693 | |
| EAN outer packing size XL | 4044941723716 | |
| EAN outer packing size XXL | - | |
| EAN outer packing size XXXL | - | |
| PZN outer packing size XS | - | |
| PZN outer packing size S | - | |
| PZN outer packing size M | - | |
| PZN outer packing size L | - | |
| PZN outer packing size XL | - | |
| PZN outer packing size XXL | - | |
| PZN outer packing size XXXL | - | |
| measures & size | | |
| length | 265 mm | |
| width | 250 mm | |
| height | 250 mm | |
| weights | | |
| size | net weight | gross weight |
| XS | 3.700 g | 4.200 g |
| S | 4.000 g | 4.500 g |
| M | 4.300 g | 4.800 g |
| L | 4.600 g | 5.100 g |
| XL | 4.900 g | 5.400 g |
| XXL | - | - |
| XXXL | - | - |

LOGISTIC DATA PALETTE

| general information | | |
|-----------------------|--------------|--------------|
| kind of palett | euro-palette | |
| measures & size | | |
| cartons per layer | 12 | |
| layers per palette | 7 | |
| height of the palette | 190 cm | |
| weights | | |
| size | net weight | gross weight |
| XS | 353 kg | 378 kg |
| S | 378 kg | 403 kg |
| M | 403 kg | 428 kg |
| L | 428 kg | 453 kg |
| XL | 454 kg | 479 kg |
| XXL | - | - |
| XXXL | - | - |